

LAMP1 Antibody

Catalog No: #32731



Package Size: #32731-1 50ul #32731-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	LAMP1 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB,IF
Species Reactivity	Human,Mouse,Rat
Specificity	The antibody detects endogenous level of total LAMP1 protein.
Immunogen Type	Recombinant Protein
Immunogen Description	Recombinant fusion protein of human LAMP1 (NP_005552.3).
Target Name	LAMP1
Other Names	LAMP1;CD107a;LAMPa;LGP120
Accession No.	Uniprot:P11279GeneID:3916
Uniprot	P11279
GeneID	3916
SDS-PAGE MW	90-120kDa
Concentration	1.0mg/ml
Formulation	PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.


Application Details

WB 1:500 - 1:2000IF 1:50 - 1:200

Images



Western blot analysis of extracts of various cell lines, using LAMP1 antibody.



Immunofluorescence analysis of NIH-3T3 cells using LAMP1 Rabbit pAb.

Background

The protein encoded by this gene is a member of a family of membrane glycoproteins. This glycoprotein provides selectins with carbohydrate ligands. It may also play a role in tumor cell metastasis.

Note: This product is for in vitro research use only