

## SMARCAD1 Antibody

Catalog No: #33088

Package Size: #33088-1 50ul #33088-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)


## Description

Product Name	SMARCAD1 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB,IF
Species Reactivity	Human,Mouse,Rat
Specificity	The antibody detects endogenous level of total SMARCAD1 protein.
Immunogen Type	Recombinant Protein
Immunogen Description	Recombinant fusion protein of human SMARCAD1 (NP_001121901.1).
Target Name	SMARCAD1
Other Names	SMARCAD1;ADERM;BASNS;ETL1;HEL1
Accession No.	Uniprot:Q9H4L7GeneID:56916
Uniprot	Q9H4L7
GeneID	56916
SDS-PAGE MW	Refer to figures
Concentration	1.0mg/ml
Formulation	PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.

## Application Details

WB□1:500 - 1:2000IF□1:50 - 1:100

## Images



Immunofluorescence analysis of A549 cells using SMARCAD1 antibody.

## Background

---

This gene encodes a member of the SNF subfamily of helicase proteins. The encoded protein plays a critical role in the restoration of heterochromatin organization and propagation of epigenetic patterns following DNA replication by mediating histone H3/H4 deacetylation. Mutations in this gene are associated with adermatoglyphia. Alternatively spliced transcript variants encoding multiple isoforms have been observed for this gene.

---

Note: This product is for in vitro research use only