

ABT1 Antibody

Catalog No: #33104



Package Size: #33104-1 50ul #33104-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	ABT1 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antibodies were purified by affinity purification using immunogen.
Applications	WB
Species Reactivity	Human,Mouse,Rat
Specificity	The antibody detects endogenous level of total ABT1 protein.
Immunogen Type	Recombinant Protein
Immunogen Description	Recombinant protein of human ABT1.
Target Name	ABT1
Other Names	hABT1;
Accession No.	Swiss-Prot:Q9ULW3NCBI Gene ID:29777
Uniprot	Q9ULW3
GeneID	29777;
SDS-PAGE MW	31KD
Concentration	1.0mg/ml
Formulation	Supplied at 1.0mg/mL in phosphate buffered saline (without Mg2+ and Ca2+), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

Application Details

WB 1:500 - 1:2000

Images



Western blot analysis of extracts of various cell lines, using ABT1 at 1:1000 dilution.

Background

Basal transcription of genes by RNA polymerase II requires the interaction of TATA-binding protein (TBP) with the core region of class II promoters. Studies in mouse suggest that the protein encoded by this gene likely activates basal transcription from class II promoters by interaction with TBP and the class II promoter DNA.

Note: This product is for in vitro research use only