

## Tau (Ab-396) Antibody

Catalog No: #33207

Package Size: #33207-1 50ul #33207-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)


## Description

Product Name	Tau (Ab-396) Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB
Species Reactivity	Hu
Specificity	The antibody detects endogenous levels of Tau protein around Serine 396.
Immunogen Type	Peptide
Immunogen Description	Synthesized non-phosphopeptide derived from human Tau around the phosphorylation site of Serine 396.
Target Name	Tau
Other Names	MAPT; Microtubule-associated protein tau; MTBT1; Neurofibrillary tangle protein; Paired helical filament-tau
Accession No.	Swiss-Prot: P10636NCBI Gene ID: 4137
Uniprot	P10636
GeneID	4137;
SDS-PAGE MW	78kd
Concentration	1.0mg/ml
Formulation	Rabbit IgG in phosphate buffered saline (without Mg <sup>2+</sup> and Ca <sup>2+</sup> ), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C


## Application Details

Western blotting: 1:500~1:3000

## Images



Western blot analysis of extracts from HeLa cells, using Tau (Ab-396) antibody #33207.



Western blot analysis of extracts from 293 cells (Lane 2), using Tau (Ab-396) antibody #33207. The lane on the left is treated with synthesized peptide.

## Background

Promotes microtubule assembly and stability, and might be involved in the establishment and maintenance of neuronal polarity. The C-terminus binds axonal microtubules while the N-terminus binds neural plasma membrane components, suggesting that tau functions as a linker protein between both. Axonal polarity is predetermined by TAU/MAPT localization (in the neuronal cell) in the domain of the cell body defined by the centrosome. The short isoforms allow plasticity of the cytoskeleton whereas the longer isoforms may preferentially play a role in its stabilization.

Puig B, et al. (2005) *Acta Neuropathol (Berl)*. 110(3):261-268.

DeGiorgis JA, et al. (2005) *J Proteome Res*. 4(2): 306-315.

Alonso Adel C, et al. (2004) *J Biol Chem*. 279(33): 34873-34881.

Note: This product is for in vitro research use only