

PEK/PERK (Ab-981) Antibody

Catalog No: #33247



Package Size: #33247-1 50ul #33247-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	PEK/PERK (Ab-981) Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB IHC IF/ICC ELISA
Species Reactivity	Hu Ms Rt
Specificity	The antibody detects endogenous levels of total PEK/PERK protein.
Immunogen Type	Peptide
Immunogen Description	Synthesized non-phosphopeptide derived from human PEK/PERK around the phosphorylation site of threonine 981 (R-H-T(p)-G-Q).
Target Name	PEK/PERK
Other Names	E2AK3; EC 2.7.11.1; EIF2AK3; Eukaryotic translation initiation factor 2-alpha kinase 3 precursor; HsPEK
Accession No.	Swiss-Prot: Q9NZJ5NCBI Gene ID: 9451
Uniprot	Q9NZJ5
GeneID	9451;
SDS-PAGE MW	125kd
Concentration	1.0mg/ml
Formulation	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

Application Details


Western blotting: 1:500~1:3000

Immunohistochemistry: 1:50~1:100


IF/ICC 1:100-1:500

ELISA 1:20000-1:40000


Images



Immunohistochemistry analysis of paraffin-embedded human colon carcinoma tissue using PEK/PERK (Ab-981) antibody #33247.



Western blot analysis of extracts from MCF-7 cells, using PEK/PERK (Ab-981) antibody #33247.



Western blot analysis of extracts from JK cells (Lane 2), using PEK/PERK (Ab-981) antibody #33247. The lane on the left is treated with synthesized peptide.

Background

Phosphorylates the alpha subunit of eukaryotic translation-initiation factor 2 (EIF2), leading to its inactivation and thus to a rapid reduction of translational initiation and repression of global protein synthesis. Serves as a critical effector of unfolded protein response (UPR)-induced G1 growth arrest due to the loss of cyclin-D1 (CCND1) By similarity.

Shi Y., J. Biol. Chem. 274:5723-5730(1999).

Sood R., Biochem. J. 346:281-293(2000).

Delepine M., Nat. Genet. 25:406-409(2000).

Note: This product is for in vitro research use only