

OVOL1 Antibody

Catalog No: #33567

Package Size: #33567-1 50ul #33567-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	OVOL1 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB
Species Reactivity	Hu Ms
Specificity	The antibody detects endogenous levels of total OVOL1 protein.
Immunogen Type	Peptide
Immunogen Description	Synthesized peptide derived from N-terminal of human OVOL1.
Target Name	OVOL1
Other Names	Putative transcription factor Ovo-like 1; hOvo1; OVOL1;
Accession No.	Swiss-Prot: O14753NCBI Gene ID: 5017
Uniprot	O14753
GeneID	5017;
SDS-PAGE MW	33kd
Concentration	1.0mg/ml
Formulation	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

Application Details

Western blotting: 1:500~1:3000

Images



Western blot analysis of extracts from K562 cells and A549 cells, using OVOL1 antibody #33567.

Background

Putative transcription factor. Involved in hair formation and spermatogenesis. May function in the differentiation and/or maintenance of the urogenital system By similarity.

The German cDNA consortium, Submitted (SEP-2004) to the EMBL/GenBank/DDBJ databases.

The MGC Project Team, Genome Res. 14:2121-2127(2004).

Chidambaram A., Mamm. Genome 8:950-951(1997).

Note: This product is for in vitro research use only