

TAF5L Antibody

Catalog No: #33615

Package Size: #33615-1 50ul #33615-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	TAF5L Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB
Species Reactivity	Hu
Specificity	The antibody detects endogenous levels of total TAF5L protein.
Immunogen Type	Peptide
Immunogen Description	Synthesized peptide derived from internal of human TAF5L.
Target Name	TAF5L
Other Names	TAF5-like RNA polymerase II p300/CBP-associated factor-associated factor 65 kDa subunit 5L; PCAF-associated factor 65 beta; PAF65-beta; TAF5L; PAF65B
Accession No.	Swiss-Prot: O75529NCBI Gene ID: 27097
Uniprot	O75529
GeneID	27097;
SDS-PAGE MW	66kd
Concentration	1.0mg/ml
Formulation	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

Application Details

Western blotting: 1:500~1:3000

Images



Western blot analysis of extracts from HT-29 cells and COS-7 cells, using TAF5L antibody #33615.

Background

Functions as a component of the PCAF complex. The PCAF complex is capable of efficiently acetylating histones in a nucleosomal context. The PCAF complex could be considered as the human version of the yeast SAGA complex.

Ogryzko V.V., Cell 94:35-44(1998).

Stanchi F., Yeast 18:69-80(2001).

Struhl K., Cell 94:1-4(1998).

Note: This product is for in vitro research use only