

ZNF682 Antibody

Catalog No: #33637

Package Size: #33637-1 50ul #33637-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	ZNF682 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB
Species Reactivity	Hu
Specificity	The antibody detects endogenous levels of total ZNF682 protein.
Immunogen Type	Peptide
Immunogen Description	Synthesized peptide derived from internal of human ZNF682.
Target Name	ZNF682
Other Names	Zinc finger protein 682; ZNF682;
Accession No.	Swiss-Prot: O95780NCBI Gene ID: 91120
Uniprot	O95780
GeneID	91120;
SDS-PAGE MW	58kd
Concentration	1.0mg/ml
Formulation	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

Application Details

Western blotting: 1:500~1:3000

Images



Western blot analysis of extracts from COLO205 cells and HUVEC cells, using ZNF682 antibody #33637.

Background

May be involved in transcriptional regulation.

Ota T., Nat. Genet. 36:40-45(2004).

Grimwood J., Nature 428:529-535(2004).

Note: This product is for in vitro research use only