

HKR1 Antibody

Catalog No: #33666

Package Size: #33666-1 50ul #33666-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description


Product Name	HKR1 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB IF
Species Reactivity	Hu
Specificity	The antibody detects endogenous levels of total HKR1 protein.
Immunogen Type	Peptide
Immunogen Description	Synthesized peptide derived from internal of human HKR1.
Target Name	HKR1
Other Names	Krueppel-related zinc finger protein 1; Protein HKR1; HKR1;
Accession No.	Swiss-Prot: P10072NCBI Gene ID: 284459
Uniprot	P10072
GeneID	284459;
SDS-PAGE MW	73kd
Concentration	1.0mg/ml
Formulation	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

Application Details

Western blotting: 1:500~1:3000


Immunofluorescence: 1:100~1:500

Images



Western blot analysis of extracts from HeLa cells, Jurkat cells and COLO205 cells, using HKR1 antibody #33666.

Immunofluorescence analysis of HUVEC cells, using HKR1 antibody #33666.



Background

May be involved in transcriptional regulation.

Oguri T., Gene 222:61-67(1998).

The MGC Project Team, Genome Res. 14:2121-2127(2004).

Ruppert J.M., Mol. Cell. Biol. 8:3104-3113(1988).

Note: This product is for in vitro research use only