

GLU2B Antibody

Catalog No: #33679



Package Size: #33679-1 50ul #33679-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description


Product Name	GLU2B Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB IF
Species Reactivity	Hu Ms
Specificity	The antibody detects endogenous levels of total GLU2B protein.
Immunogen Type	Peptide
Immunogen Description	Synthesized peptide derived from internal of human GLU2B.
Target Name	GLU2B
Other Names	80K-H; AGE-R2; AGE-binding receptor 2; G19P1; GLU2B
Accession No.	Swiss-Prot: P14314NCBI Gene ID: 5589
Uniprot	P14314
GeneID	5589;
SDS-PAGE MW	59kd
Concentration	1.0mg/ml
Formulation	Rabbit IgG in phosphate buffered saline (without Mg2+ and Ca2+), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

Application Details

Western blotting: 1:500~1:3000


Immunofluorescence: 1:100~1:500

Images



Western blot analysis of extracts from HepG2, 293 and HUVEC cells, using GLU2B antibody #33679.

Immunofluorescence analysis of HeLa cells, using GLU2B antibody #33679.



Background

Regulatory subunit of glucosidase II.

Ophoff R.A., Eur. J. Hum. Genet. 4:321-328(1996).

Pelletier M.F., Glycobiology 10:815-827(2000).

Gevaert K., Nat. Biotechnol. 21:566-569(2003).

Note: This product is for in vitro research use only