

MYF6 Antibody

Catalog No: #33719



Package Size: #33719-1 50ul #33719-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	MYF6 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB
Species Reactivity	Hu Ms Rt
Specificity	The antibody detects endogenous levels of total MYF6 protein.
Immunogen Type	Peptide
Immunogen Description	Synthesized peptide derived from internal of human MYF6.
Target Name	MYF6
Other Names	MRF4; Myf6; myogenic factor 6;
Accession No.	Swiss-Prot: P23409NCBI Gene ID: 4618
Uniprot	P23409
GeneID	4618;
SDS-PAGE MW	26kd
Concentration	1.0mg/ml
Formulation	Rabbit IgG in phosphate buffered saline (without Mg2+ and Ca2+), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

Application Details

Western blotting: 1:500~1:3000

Images



Western blot analysis of extracts from 293 cells and HepG2 cells, using MYF6 antibody #33719.

Background

Involved in muscle differentiation (myogenic factor). Induces fibroblasts to differentiate into myoblasts. Probable sequence specific DNA-binding protein.

Braun T., EMBO J. 9:821-831(1990).

The MGC Project Team, Genome Res. 14:2121-2127(2004).

Kerst B., Hum. Mutat. 12:220-220(1998).

Note: This product is for in vitro research use only