

HOXD12 Antibody

Catalog No: #33758

Package Size: #33758-1 50ul #33758-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	HOXD12 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB
Species Reactivity	Hu Ms
Specificity	The antibody detects endogenous levels of total HOXD12 protein.
Immunogen Type	Peptide
Immunogen Description	Synthesized peptide derived from internal of human HOXD12.
Target Name	HOXD12
Other Names	Homeobox protein Hox-D12; Hox-4H; HOXD12; HOX4H;
Accession No.	Swiss-Prot: P35452NCBI Gene ID: 3238
Uniprot	P35452
GeneID	3238;
SDS-PAGE MW	28kd
Concentration	1.0mg/ml
Formulation	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

Application Details

Western blotting: 1:500~1:3000

Images



Western blot analysis of extracts from K562, JK and 3T3 cells using HOXD12 antibody #33758.

Background

Sequence-specific transcription factor which is part of a developmental regulatory system that provides cells with specific positional identities on the anterior-posterior axis.

Miwa N., Submitted (MAY-1999) to the EMBL/GenBank/DDBJ databases.

D'Esposito M., Genomics 10:43-50(1991).

Note: This product is for in vitro research use only