

ETV4 Antibody

Catalog No: #33772

Package Size: #33772-1 50ul #33772-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	ETV4 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	IHC
Species Reactivity	Hu Ms
Specificity	The antibody detects endogenous levels of total ETV4 protein.
Immunogen Type	Peptide
Immunogen Description	Synthesized peptide derived from internal of human ETV4.
Target Name	ETV4
Other Names	ETS translocation variant 4; Adenovirus E1A enhancer-binding protein; E1A-F; ETV4; E1AF
Accession No.	Swiss-Prot: P43268NCBI Gene ID: 2118
Uniprot	P43268
GeneID	2118;
SDS-PAGE MW	54kd
Concentration	1.0mg/ml
Formulation	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

Application Details

Immunohistochemistry: 1:50~1:100

Images



Immunohistochemistry analysis of paraffin-embedded human breast carcinoma tissue, using ETV4 antibody #33772.

Background

Transcriptional activator that binds to the enhancer of the adenovirus E1A gene; the core-binding sequence is 5'[AC]GGA[AT]GT-3'.

Friedman L.S., Genomics 25:256-263(1995).

Coutte L., Gene 240:201-207(1999).

Higashino F., J. Biol. Chem. 268:547-553(1993).

Note: This product is for in vitro research use only