

SCAND1 Antibody

Catalog No: #33828



Package Size: #33828-1 50ul #33828-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description


Product Name	SCAND1 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB IF
Species Reactivity	Hu
Specificity	The antibody detects endogenous levels of total SCAND1 protein.
Immunogen Type	Peptide
Immunogen Description	Synthesized peptide derived from internal of human SCAND1.
Target Name	SCAND1
Other Names	SCAN domain-containing protein 1; SCAND1; SDP1;
Accession No.	Swiss-Prot: P57086NCBI Gene ID: 51282
Uniprot	P57086
GeneID	51282;
SDS-PAGE MW	70,25kd
Concentration	1.0mg/ml
Formulation	Rabbit IgG in phosphate buffered saline (without Mg2+ and Ca2+), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

Application Details

Western blotting: 1:500~1:3000


Immunofluorescence: 1:100~1:500

Images



Western blot analysis of extracts from Jurkat cells, using SCAND1 antibody #33828.

Immunofluorescence analysis of HepG2 cells, using SCAND1 antibody #33828.



Background

May regulate transcriptional activity.

Schumacher C., J. Biol. Chem. 275:17173-17179(2000).

Deloukas P., Nature 414:865-871(2001).

The MGC Project Team, Genome Res. 14:2121-2127(2004).

Note: This product is for in vitro research use only