

SFRS17A Antibody

Catalog No: #33844

Package Size: #33844-1 50ul #33844-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	SFRS17A Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB
Species Reactivity	Hu Ms
Specificity	The antibody detects endogenous levels of total SFRS17A protein.
Immunogen Type	Peptide
Immunogen Description	Synthesized peptide derived from C-terminal of human SFRS17A.
Target Name	SFRS17A
Other Names	A-kinase anchor protein 17A; 721P; B-lymphocyte antigen; Protein XE7; Protein kinase A-anchoring protein 17A
Accession No.	Swiss-Prot: Q02040NCBI Gene ID: 8227
Uniprot	Q02040
GeneID	8227;
SDS-PAGE MW	80kd
Concentration	1.0mg/ml
Formulation	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

Application Details

Western blotting: 1:500~1:3000

Images



Western blot analysis of extracts from NIH-3T3 cells, using SFRS17A antibody #33844.

Background

Splice factor regulating alternative splice site selection for certain mRNA precursors. Mediates regulation of pre-mRNA splicing in a PKA-dependent manner.

Ellison J.W., Hum. Mol. Genet. 1:691-696(1992).

Voland J.R., Proc. Natl. Acad. Sci. U.S.A. 89:10425-10429(1992).

Ross M.T., Nature 434:325-337(2005).

Note: This product is for in vitro research use only