

MAP3K10 Antibody

Catalog No: #33848

Package Size: #33848-1 50ul #33848-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	MAP3K10 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	IHC
Species Reactivity	Hu
Specificity	The antibody detects endogenous levels of total MAP3K10 protein.
Immunogen Type	Peptide
Immunogen Description	Synthesized peptide derived from internal of human MAP3K10.
Target Name	MAP3K10
Other Names	EC 2.7.11.25; Mixed lineage kinase 2; Protein kinase MST;
Accession No.	Swiss-Prot: Q02779NCBI Gene ID: 4294
Uniprot	Q02779
GeneID	4294;
SDS-PAGE MW	104kd
Concentration	1.0mg/ml
Formulation	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

Application Details

Immunohistochemistry: 1:50~1:100

Images



Immunohistochemistry analysis of paraffin-embedded human brain tissue, using MAP3K10 antibody #33848.

Background

Activates the JUN N-terminal pathway By similarity.

Dorow D.S., Eur. J. Biochem. 234:492-500(1995).

Katoh M., Oncogene 10:1447-1451(1995).

Dorow D.S., Eur. J. Biochem. 213:701-710(1993).

Note: This product is for in vitro research use only