

CEBPZ Antibody

Catalog No: #33853

Package Size: #33853-1 50ul #33853-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	CEBPZ Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB
Species Reactivity	Hu
Specificity	The antibody detects endogenous levels of total CEBPZ protein.
Immunogen Type	Peptide
Immunogen Description	Synthesized peptide derived from internal of human CEBPZ.
Target Name	CEBPZ
Other Names	CBF; CBF2; CCAAT-binding factor; CCAAT-box-binding transcription factor;
Accession No.	Swiss-Prot: Q03701NCBI Gene ID: 10153
Uniprot	Q03701
GeneID	10153;
SDS-PAGE MW	120kd
Concentration	1.0mg/ml
Formulation	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

Application Details

Western blotting: 1:500~1:3000

Images



Western blot analysis of extracts from COLO cells, using CEBPZ antibody #33853.

Background

Stimulates transcription from the HSP70 promoter.

Lum L., Mol. Cell. Biol. 10:6709-6717(1990).

Hillier L.W., Nature 434:724-731(2005).

Olsen J.V., Cell 127:635-648(2006).

Note: This product is for in vitro research use only