

FLI1 Antibody

Catalog No: #21433

Package Size: #21433-1 50ul #21433-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description


Product Name	FLI1 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antibodies were produced by immunizing rabbits with synthetic peptide and KLH conjugates. Antibodies were purified by affinity-chromatography using epitope-specific peptide.
Applications	WB
Species Reactivity	Hu
Specificity	The antibody detects endogenous levels of total FLI1 protein.
Immunogen Type	Peptide-KLH
Immunogen Description	Peptide sequence around aa.251~255(S-K-N-T-E)derived from Human FLI1.
Target Name	FLI1
Other Names	EWSR2; SIC-1;
Accession No.	Swiss-Prot: Q01543NCBI Protein: NP_002008.2
Uniprot	Q01543
GeneID	2313;
Concentration	1.0mg/ml
Formulation	Supplied at 1.0mg/mL in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.

Application Details

Predicted MW: 51kd

Western blotting: 1:500~1:1000

Images



Western blot analysis of extracts from 293 and HeLa cells using FLI1 Antibody #21433.

Background

Sequence-specific transcriptional activator. Recognizes the DNA sequence 5'-C[CA]GGAAGT-3'.

Delattre O., Zucman J., Plougastel B., Desmaze C. (1992) Nature 359:162-165 Watson D.K., Smyth F.E., Thompson D.M., (1992) Cell Growth Differ. 3:705-713

Note: This product is for in vitro research use only