

CNKR2 Antibody

Catalog No: #33998

Package Size: #33998-1 50ul #33998-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	CNKR2 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB
Species Reactivity	Hu Ms Rt
Specificity	The antibody detects endogenous levels of total CNKR2 protein.
Immunogen Type	Peptide
Immunogen Description	Synthesized peptide derived from internal of human CNKR2.
Target Name	CNKR2
Other Names	Connector enhancer of kinase suppressor of ras 2; Connector enhancer of KSR2; CNK2; CNKSR2; CNK2
Accession No.	Swiss-Prot: Q8WXI2NCBI Gene ID: 22866
Uniprot	Q8WXI2
GeneID	22866;
SDS-PAGE MW	117kd
Concentration	1.0mg/ml
Formulation	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

Application Details

Western blotting: 1:500~1:3000

Images



Western blot analysis of extracts from Jurkat cells, using CNKR2 antibody #33998.

Background

May function as an adapter protein or regulator of Ras signaling pathways.

Lanigan T.M., FASEB J. 17:2048-2060(2003).

Nagase T., DNA Res. 5:355-364(1998).

Ross M.T., Nature 434:325-337(2005).

Note: This product is for in vitro research use only