

SEPT1 Antibody

Catalog No: #34001

Package Size: #34001-1 50ul #34001-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	SEPT1 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB
Species Reactivity	Hu Ms Rt
Specificity	The antibody detects endogenous levels of total SEPT1 protein.
Immunogen Type	Peptide
Immunogen Description	Synthesized peptide derived from internal of human SEPT1.
Target Name	SEPT1
Other Names	Septin-1; LARP; Serologically defined breast cancer antigen NY-BR-24; Peanut-like protein 3; SEPT1
Accession No.	Swiss-Prot: Q8WYJ6NCBI Gene ID: 1731
Uniprot	Q8WYJ6
GeneID	1731;
SDS-PAGE MW	41kd
Concentration	1.0mg/ml
Formulation	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

Application Details

Western blotting: 1:500~1:3000

Images



Western blot analysis of extracts from Jurkat cells, using SEPT1 antibody #34001.

Background

Filament-forming cytoskeletal GTPase By similarity. May play a role in cytokinesis Potential.

Tu Q., Submitted (AUG-1998) to the EMBL/GenBank/DDBJ databases.

The MGC Project Team, Genome Res. 14:2121-2127(2004).

Scanlan M.J., Cancer Immun. 1:4-4(2001).

Note: This product is for in vitro research use only