

MOBKL2A/B Antibody

Catalog No: #34025

Package Size: #34025-1 50ul #34025-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	MOBKL2A/B Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB
Species Reactivity	Hu Ms
Specificity	The antibody detects endogenous levels of total MOBKL2A/B protein.
Immunogen Type	Peptide
Immunogen Description	Synthesized peptide derived from internal of human MOBKL2A/B.
Target Name	MOBKL2A/B
Other Names	Mps one binder kinase activator-like 2A; MOB-LAK; Mob1 homolog 2A; Protein Mob3A;
Accession No.	Swiss-Prot: Q96BX8/Q86TA1NCBI Gene ID: 126308/79817
Uniprot	Q96BX8
GeneID	126308;
SDS-PAGE MW	25kd
Concentration	1.0mg/ml
Formulation	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

Application Details

Western blotting: 1:500~1:3000

Images



Western blot analysis of extracts from 3T3 cells, using MOBKL2A/B antibody #34025.

Background

May regulate the activity of kinases By similarity.

Ota T., Nat. Genet. 36:40-45(2004) [PubMed: 14702039].

Grimwood J., Nature 428:529-535(2004) [PubMed: 15057824].

Note: This product is for in vitro research use only