

ZNF287 Antibody

Catalog No: #34071

Package Size: #34071-1 50ul #34071-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description


Product Name	ZNF287 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB IHC
Species Reactivity	Hu Ms
Specificity	The antibody detects endogenous levels of total ZNF287 protein.
Immunogen Type	Peptide
Immunogen Description	Synthesized peptide derived from internal of human ZNF287.
Target Name	ZNF287
Other Names	Zinc finger protein with KRAB and SCAN domains 13; Zinc finger protein 287; ZNF287; ZKSCAN13;
Accession No.	Swiss-Prot: Q9HBT7NCBI Gene ID: 57336
Uniprot	Q9HBT7
GeneID	57336;
SDS-PAGE MW	90kd
Concentration	1.0mg/ml
Formulation	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

Application Details

Western blotting: 1:500~1:3000


Immunohistochemistry: 1:50~1:100

Images



Western blot analysis of extracts from NIH-3T3 cells, using ZNF287 antibody #34071.

Immunohistochemistry analysis of paraffin-embedded human brain tissue using ZNF287 antibody #34071.



Background

May be involved in transcriptional regulation.

Romero-Pastrana F., Submitted (DEC-1999) to the EMBL/GenBank/DDBJ databases.

Ebert L., Submitted (JUN-2004) to the EMBL/GenBank/DDBJ databases.

The MGC Project Team, Genome Res. 14:2121-2127(2004).

Note: This product is for in vitro research use only