

ARMX1 Antibody

Catalog No: #34090

Package Size: #34090-1 50ul #34090-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	ARMX1 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB
Species Reactivity	Hu Rt
Specificity	The antibody detects endogenous levels of total ARMX1 protein.
Immunogen Type	Peptide
Immunogen Description	Synthesized peptide derived from internal of human ARMX1.
Target Name	ARMX1
Other Names	Armadillo repeat-containing X-linked protein 1; Protein ALEX1; ARM protein lost in epithelial cancers on chromosome X 1; ARMCX1; ALEX1
Accession No.	Swiss-Prot: Q9P291NCBI Gene ID: 51309
Uniprot	Q9P291
GeneID	51309;
SDS-PAGE MW	49kd
Concentration	1.0mg/ml
Formulation	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

Application Details

Western blotting: 1:500~1:3000

Images



Western blot analysis of extracts from rat brain cells, using ARMX1 antibody #34090.

Background

Kurochkin I.V., Biochem. Biophys. Res. Commun. 280:340-347(2001).

Ota T., Nat. Genet. 36:40-45(2004).

Ross M.T., Nature 434:325-337(2005).

Note: This product is for in vitro research use only