

ZNF337 Antibody

Catalog No: #34133

Package Size: #34133-1 50ul #34133-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	ZNF337 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB
Species Reactivity	Hu
Specificity	The antibody detects endogenous levels of total ZNF337 protein.
Immunogen Type	Peptide
Immunogen Description	Synthesized peptide derived from internal of human ZNF337.
Target Name	ZNF337
Other Names	Zinc finger protein 337; ZN337; ZNF337;
Accession No.	Swiss-Prot: Q9Y3M9NCBI Gene ID: 26152
Uniprot	Q9Y3M9
GeneID	26152;
SDS-PAGE MW	86kd
Concentration	1.0mg/ml
Formulation	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

Application Details

Western blotting: 1:500~1:3000

Images



Western blot analysis of extracts from HeLa cells, using ZNF337 antibody #34133.

Background

May be involved in transcriptional regulation.

Deloukas P., Nature 414:865-871(2001).

The MGC Project Team; Genome Res. 14:2121-2127(2004).

The German cDNA consortium; Submitted (MAY-1999) to the EMBL/GenBank/DDBJ databases.

Note: This product is for in vitro research use only