

TF3C3 Antibody

Catalog No: #34137

Package Size: #34137-1 50ul #34137-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	TF3C3 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB
Species Reactivity	Hu Ms
Specificity	The antibody detects endogenous levels of total TF3C3 protein.
Immunogen Type	Peptide
Immunogen Description	Synthesized peptide derived from internal of human TF3C3.
Target Name	TF3C3
Other Names	TF3C-gamma; TF3C3; TFIIIC 102 kDa subunit; TFIIIC102; TFIIICgamma
Accession No.	Swiss-Prot: Q9Y5Q9NCBI Gene ID: 9330
Uniprot	Q9Y5Q9
GeneID	9330;
SDS-PAGE MW	101kd
Concentration	1.0mg/ml
Formulation	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

Application Details

Western blotting: 1:500~1:3000

Images



Western blot analysis of extracts from 293 cells, using TF3C3 antibody #34137.

Background

Involved in RNA polymerase III-mediated transcription. Integral, tightly associated component of the DNA-binding TFIIIC2 subcomplex that directly binds tRNA and virus-associated RNA promoters.

Ota T., Nat. Genet. 36:40-45(2004).

The German cDNA consortium; Submitted (JAN-2000) to the EMBL/GenBank/DDBJ databases.

Hsieh Y.-J., Mol. Cell. Biol. 19:4944-4952(1999).

Note: This product is for in vitro research use only