

ARRDC2 Antibody

Catalog No: #34168

Package Size: #34168-1 50ul #34168-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	ARRDC2 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	IHC
Species Reactivity	Hu Ms
Specificity	The antibody detects endogenous levels of total ARRDC2 protein.
Immunogen Type	Peptide
Immunogen Description	Synthesized peptide derived from internal of human ARRDC2.
Target Name	ARRDC2
Other Names	Arrestin domain-containing protein 2; ARRDC2; PP2703;
Accession No.	Swiss-Prot: Q8TBH0NCBI Gene ID: 27106
Uniprot	Q8TBH0
GeneID	27106;
SDS-PAGE MW	44kd
Concentration	1.0mg/ml
Formulation	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

Application Details

Immunohistochemistry: 1:50~1:100

Images



Immunohistochemistry analysis of paraffin-embedded human brain tissue, using ARRDC2 antibody #34168.

Background

Belongs to the arrestin family. 2 isoforms of the human protein are produced by alternative splicing.

Ota T., Nat. Genet. 36:40-45(2004).

Wan D., Proc. Natl. Acad. Sci. U.S.A. 101:15724-15729(2004).

The MGC Project Team, Genome Res. 14:2121-2127(2004).

Note: This product is for in vitro research use only