

C56D2 Antibody

Catalog No: #34171

Package Size: #34171-1 50ul #34171-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	C56D2 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB
Species Reactivity	Hu
Specificity	The antibody detects endogenous levels of total C56D2 protein.
Immunogen Type	Peptide
Immunogen Description	Synthesized peptide derived from internal of human C56D2.
Target Name	C56D2
Other Names	Cytochrome b561 domain-containing protein 2; Putative tumor suppressor protein 101F6; CYB561D2; 101F6; LUCA12.2
Accession No.	Swiss-Prot: O14569NCBI Gene ID: 11068
Uniprot	O14569
GeneID	11068;
SDS-PAGE MW	24kd
Concentration	1.0mg/ml
Formulation	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

Application Details

Western blotting: 1:500~1:3000

Images



Western blot analysis of extracts from Jurkat cells, using C56D2 antibody #34171.

Background

Bader S., Submitted (JAN-1998) to the EMBL/GenBank/DDBJ databases.

Muzny D.M., Nature 440:1194-1198(2006).

The MGC Project Team; Genome Res. 14:2121-2127(2004).

Note: This product is for in vitro research use only