

DNAJB4 Antibody

Catalog No: #34246

Package Size: #34246-1 50ul #34246-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	DNAJB4 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB
Species Reactivity	Hu Ms Rt
Specificity	The antibody detects endogenous levels of total DNAJB4 protein.
Immunogen Type	Peptide
Immunogen Description	Synthesized peptide derived from internal of human DNAJB4.
Target Name	DNAJB4
Other Names	DJB4; DnaJ homolog subfamily B member 4; DNAJW; DNJB4; Heat shock 40 kDa protein 1 homolog
Accession No.	Swiss-Prot: Q9UDY4NCBI Gene ID: 11080
Uniprot	Q9UDY4
GeneID	11080;
SDS-PAGE MW	38kd
Concentration	1.0mg/ml
Formulation	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

Application Details

Western blotting: 1:500~1:3000

Images



Western blot analysis of extracts from HepG2 cells, using DNAJB4 antibody #34246.

Background

Probable chaperone.

Hoe K.L., Biochim. Biophys. Acta 1383:4-8(1998).

The MGC Project Team; Genome Res. 14:2121-2127(2004).

Rush J., Nat. Biotechnol. 23:94-101(2005).

Note: This product is for in vitro research use only