

IL20RB Antibody

Catalog No: #34254

Package Size: #34254-1 50ul #34254-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	IL20RB Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB
Species Reactivity	Hu Ms
Specificity	The antibody detects endogenous levels of total IL20RB protein.
Immunogen Type	Peptide
Immunogen Description	Synthesized peptide derived from internal of human IL20RB.
Target Name	IL20RB
Other Names	DIRS1; FNDC6; I20RB; IL-20R-beta; IL-20R2
Accession No.	Swiss-Prot: Q6UXL0NCBI Gene ID: 53833
Uniprot	Q6UXL0
GeneID	53833;
SDS-PAGE MW	35kd
Concentration	1.0mg/ml
Formulation	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

Application Details

Western blotting: 1:500~1:3000

Images



Western blot analysis of extracts from HeLa cells, using IL20RB antibody #34254.

Background

The IL20RA/IL20RB dimer is a receptor for IL19, IL20 and IL24. The IL22RA1/IL20RB dimer is a receptor for IL20 and IL24.

Clark H.F., Genome Res. 13:2265-2270(2003).

The MGC Project Team; Genome Res. 14:2121-2127(2004).

Zhang Z., Protein Sci. 13:2819-2824(2004).

Note: This product is for in vitro research use only