

TPIP1 Antibody

Catalog No: #34300

Package Size: #34300-1 50ul #34300-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description


Product Name	TPIP1 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	IHC IF
Species Reactivity	Hu
Specificity	The antibody detects endogenous levels of total TPIP1 protein.
Immunogen Type	Peptide
Immunogen Description	Synthesized peptide derived from internal of human TPIP1.
Target Name	TPIP1
Other Names	p53-regulated apoptosis-inducing protein 1; P53aip1; TP53AIP1;
Accession No.	Swiss-Prot: Q9HCN2NCBI Gene ID: 63970
Uniprot	Q9HCN2
GeneID	63970;
SDS-PAGE MW	13kd
Concentration	1.0mg/ml
Formulation	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

Application Details

Immunohistochemistry: 1:50~1:100


Immunofluorescence: 1:100~1:500

Images



Immunohistochemistry analysis of paraffin-embedded human thymus gland tissue using TPIP1 antibody #34300.

Immunofluorescence analysis of A549 cells, using TPIP1 antibody #34300.



Background

May play an important role in mediating p53/TP53-dependent apoptosis.

Oda K., Cell 102:849-862(2000).

Zhou G., Submitted (MAY-2003) to the EMBL/GenBank/DDBJ databases.

The MGC Project Team, Genome Res. 14:2121-2127(2004).

Note: This product is for in vitro research use only