

MP68 Antibody

Catalog No: #34343

Package Size: #34343-1 50ul #34343-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	MP68 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	IHC
Species Reactivity	Hu Ms
Specificity	The antibody detects endogenous levels of total MP68 protein.
Immunogen Type	Peptide
Immunogen Description	Synthesized peptide derived from internal of human MP68.
Target Name	MP68
Other Names	6.8 kDa mitochondrial proteolipid;
Accession No.	Swiss-Prot: P56378NCBI Gene ID: 9556
Uniprot	P56378
GeneID	9556;
SDS-PAGE MW	6kd
Concentration	1.0mg/ml
Formulation	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

Application Details

Immunohistochemistry: 1:50~1:100

Images



Immunohistochemistry analysis of paraffin-embedded human prostate carcinoma tissue using MP68 antibody #34343.

Background

Belongs to the small mitochondrial proteolipid family.

Li W.B., Submitted (JAN-2003) to the EMBL/GenBank/DDBJ databases.

Mao M., Proc. Natl. Acad. Sci. U.S.A. 95:8175-8180(1998).

Note: This product is for in vitro research use only