

RPLP2 Antibody

Catalog No: #34344

Package Size: #34344-1 50ul #34344-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	RPLP2 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	IHC
Species Reactivity	Hu Ms
Specificity	The antibody detects endogenous levels of total RPLP2 protein.
Immunogen Type	Peptide
Immunogen Description	Synthesized peptide derived from internal of human RPLP2.
Target Name	RPLP2
Other Names	60S acidic ribosomal protein P2; Renal carcinoma antigen NY-REN-44;
Accession No.	Swiss-Prot: P05387NCBI Gene ID: 6181
Uniprot	P05387
GeneID	6181;
SDS-PAGE MW	12kd
Concentration	1.0mg/ml
Formulation	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

Application Details

Immunohistochemistry: 1:50~1:100

Images



Immunohistochemistry analysis of paraffin-embedded human prostate carcinoma tissue using RPLP2 antibody #34344.

Background

Plays an important role in the elongation step of protein synthesis. HAMAP-Rule MF_01478

Yoshihama M., Genome Res. 12:379-390(2002).

Giorgianni F., Proteomics 4:587-598(2004).

Kim J.-E., J. Proteome Res. 4:1339-1346(2005).

Note: This product is for in vitro research use only