

RPL10L Antibody

Catalog No: #34345

Package Size: #34345-1 50ul #34345-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	RPL10L Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB
Species Reactivity	Hu
Specificity	The antibody detects endogenous levels of total RPL10L protein.
Immunogen Type	Peptide
Immunogen Description	Synthesized peptide derived from internal of human RPL10L.
Target Name	RPL10L
Other Names	60S ribosomal protein L10-like; RL10L;
Accession No.	Swiss-Prot: Q96L21NCBI Gene ID: 140801
Uniprot	Q96L21
GeneID	140801;
SDS-PAGE MW	20kd
Concentration	1.0mg/ml
Formulation	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

Application Details

Western blotting: 1:500~1:3000

Images



Western blot analysis of extracts from HeLa cells, Jurkat cells and COS7 cells, using RPL10L antibody #34345.

Background

May play a role in compensating for the inactivated X-linked gene during spermatogenesis.

Uechi T., Nucleic Acids Res. 30:5369-5375(2002).

The MGC Project Team; Genome Res. 14:2121-2127(2004).

Note: This product is for in vitro research use only