

RPL12 Antibody

Catalog No: #34346

Package Size: #34346-1 50ul #34346-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	RPL12 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB
Species Reactivity	Hu Ms
Specificity	The antibody detects endogenous levels of total RPL12 protein.
Immunogen Type	Peptide
Immunogen Description	Synthesized peptide derived from internal of human RPL12.
Target Name	RPL12
Other Names	RL12; RPL12;
Accession No.	Swiss-Prot: P30050NCBI Gene ID: 6136
Uniprot	P30050
GeneID	6136;
SDS-PAGE MW	18kd
Concentration	1.0mg/ml
Formulation	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

Application Details

Western blotting: 1:500~1:3000

Images



Western blot analysis of extracts from HepG2 cells and HeLa cells, using RPL12 antibody #34346.

Background

Binds directly to 26S ribosomal RNA By similarity. HAMAP-Rule MF_00736

Chu W., Nucleic Acids Res. 21:749-749(1993).

Humphray S.J., Nature 429:369-374(2004).

The MGC Project Team; Genome Res. 14:2121-2127(2004).

Note: This product is for in vitro research use only