

RPL3 Antibody

Catalog No: #34353

Package Size: #34353-1 50ul #34353-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description


Product Name	RPL3 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB IHC
Species Reactivity	Hu Ms Rt
Specificity	The antibody detects endogenous levels of total RPL3 protein.
Immunogen Type	Peptide
Immunogen Description	Synthesized peptide derived from internal of human RPL3.
Target Name	RPL3
Other Names	60S ribosomal protein L3; HIV-1 TAR RNA binding protein B; RL3; RPL3; TARBP-B
Accession No.	Swiss-Prot: P39023NCBI Gene ID: 6122
Uniprot	P39023
GeneID	6122;
SDS-PAGE MW	46kd
Concentration	1.0mg/ml
Formulation	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

Application Details

Western blotting: 1:500~1:3000


Immunohistochemistry: 1:50~1:100

Images



Western blot analysis of extracts from HT-29 cells, using RPL3 antibody #34353.

Immunohistochemistry analysis of paraffin-embedded human breast carcinoma tissue using RPL3 antibody #34353.



Background

The L3 protein is a component of the large subunit of cytoplasmic ribosomes.

Leffers H., Submitted (JUN-1993) to the EMBL/GenBank/DDBJ databases.

Shichijo S., Submitted (MAY-2001) to the EMBL/GenBank/DDBJ databases.

Collins J.E., Genome Biol. 5:RESEARCH84.1-RESEARCH84.11(2004).

Note: This product is for in vitro research use only