

ABHD14A Antibody

Catalog No: #34368

Package Size: #34368-1 50ul #34368-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	ABHD14A Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB
Species Reactivity	Hu Ms
Specificity	The antibody detects endogenous levels of total ABHD14A protein.
Immunogen Type	Peptide
Immunogen Description	Synthesized peptide derived from internal of human ABHD14A.
Target Name	ABHD14A
Other Names	Abhydrolase domain-containing protein 14A; EC 3.-.-.-;
Accession No.	Swiss-Prot: Q9BUJ0NCBI Gene ID: 25864
Uniprot	Q9BUJ0
GeneID	25864;
SDS-PAGE MW	30kd
Concentration	1.0mg/ml
Formulation	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

Application Details

Western blotting: 1:500~1:3000

Images



Western blot analysis of extracts from RAW264.7 cells, using ABHD14A antibody #34368.

Background

Possible role in granule neuron development By similarity.

Clark H.F., Genome Res. 13:2265-2270(2003).

The German cDNA consortium; Submitted (SEP-2004) to the EMBL/GenBank/DDBJ databases.

Genome Res. 14:2121-2127(2004).

Note: This product is for in vitro research use only