

ABHD6 Antibody

Catalog No: #34371

Package Size: #34371-1 50ul #34371-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description


Product Name	ABHD6 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB IF
Species Reactivity	Hu Ms Rt
Specificity	The antibody detects endogenous levels of total ABHD6 protein.
Immunogen Type	Peptide
Immunogen Description	Synthesized peptide derived from internal of human ABHD6.
Target Name	ABHD6
Other Names	Abhydrolase domain-containing protein 6; EC 3.-.-.-; ABHD6;
Accession No.	Swiss-Prot: Q9BV23NCBI Gene ID: 57406
Uniprot	Q9BV23
GeneID	57406;
SDS-PAGE MW	38kd
Concentration	1.0mg/ml
Formulation	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

Application Details

Western blotting: 1:500~1:3000


Immunofluorescence: 1:100~1:500

Images



Western blot analysis of extracts from COLO cells, using ABHD6 antibody #34371.

Immunofluorescence analysis of A549 cells, using ABHD6 antibody #34371.



Background

Has 2-arachidonoylglycerol hydrolase activity By similarity. May be a regulator of endocannabinoid signaling pathways By similarity.

Ota T., Nat. Genet. 36:40-45(2004).

The MGC Project Team; Genome Res. 14:2121-2127(2004).

Note: This product is for in vitro research use only