

ACAD10 Antibody

Catalog No: #34377

Package Size: #34377-1 50ul #34377-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	ACAD10 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	IHC
Species Reactivity	Hu
Specificity	The antibody detects endogenous levels of total ACAD10 protein.
Immunogen Type	Peptide
Immunogen Description	Synthesized peptide derived from internal of human ACAD10.
Target Name	ACAD10
Other Names	Acyl-CoA dehydrogenase family member; 10ACAD-10;
Accession No.	Swiss-Prot: Q6JQN1NCBI Gene ID: 80724
Uniprot	Q6JQN1
GeneID	80724;
SDS-PAGE MW	119kd
Concentration	1.0mg/ml
Formulation	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

Application Details

Immunohistochemistry: 1:50~1:100

Images



Immunohistochemistry analysis of paraffin-embedded human heart tissue using ACAD10 antibody #34377.

Background

Acyl-CoA dehydrogenase only active with R- and S-2-methyl-C15-CoA.

Ye X., Mol. Biol. Rep. 31:191-195(2004).

Ota T., Nat. Genet. 36:40-45(2004).

The MGC Project Team; Genome Res. 14:2121-2127(2004).

Note: This product is for in vitro research use only