

ADAM32 Antibody

Catalog No: #34386

Package Size: #34386-1 50ul #34386-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	ADAM32 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB
Species Reactivity	Hu
Specificity	The antibody detects endogenous levels of total ADAM32 protein.
Immunogen Type	Peptide
Immunogen Description	Synthesized peptide derived from Internal of human ADAM32.
Target Name	ADAM32
Other Names	a disintegrin and metalloprotease domain 32; a disintegrin and metalloproteinase domain 32; ADA32; ADAM 32; ADAM metalloproteinase domain 32
Accession No.	Swiss-Prot: C14296NCBI Gene ID: 203102
Uniprot	Q8TC27
GeneID	203102;
SDS-PAGE MW	88kd
Concentration	1.0mg/ml
Formulation	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

Application Details

Western blotting: 1:500~1:3000

Images



Western blot analysis of extracts from 293 cells, using ADAM32 antibody #34386.

Background

May play a role in sperm development and fertilization This is a non-catalytic metalloprotease-like protein.

Clark H.F., Genome Res. 13:2265-2270(2003).

Nusbaum C., Nature 439:331-335(2006).

Zhang Z., Protein Sci. 13:2819-2824(2004).

Note: This product is for in vitro research use only