

AKR1CL1 Antibody

Catalog No: #34401

Package Size: #34401-1 50ul #34401-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	AKR1CL1 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	IF
Species Reactivity	Hu
Specificity	The antibody detects endogenous levels of total AKR1CL1 protein.
Immunogen Type	Peptide
Immunogen Description	Synthesized peptide derived from internal of human AKR1CL1.
Target Name	AKR1CL1
Other Names	Aldo-keto reductase family 1 member C-like protein 1; EC 1.1.1.-;
Accession No.	Swiss-Prot: Q5T2L2NCBI Gene ID: 340811
Uniprot	Q5T2L2
GeneID	340811;
SDS-PAGE MW	35kd
Concentration	1.0mg/ml
Formulation	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

Application Details

Immunofluorescence: 1:100~1:500

Images



Immunofluorescence analysis of A549 cells, using AKR1CL1 antibody #34401.

Background

Belongs to the aldo/keto reductase family.

Ota T., Nat. Genet. 36:40-45(2004).

Deloukas P., Nature 429:375-381(2004).

The MGC Project Team;

Genome Res. 14:2121-2127(2004).

Note: This product is for in vitro research use only