

APBB2 Antibody

Catalog No: #34412

Package Size: #34412-1 50ul #34412-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description


Product Name	APBB2 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB IHC
Species Reactivity	Hu Ms
Specificity	The antibody detects endogenous levels of total APBB2 protein.
Immunogen Type	Peptide
Immunogen Description	Synthesized peptide derived from internal of human APBB2.
Target Name	APBB2
Other Names	amyloid beta A4 precursor protein-binding family B member 2; Fe65-like protein;
Accession No.	Swiss-Prot: Q92870NCBI Gene ID: 323
Uniprot	Q92870
GeneID	323;
SDS-PAGE MW	83kd
Concentration	1.0mg/ml
Formulation	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

Application Details

Western blotting: 1:500~1:3000


Immunohistochemistry: 1:50~1:100

Images



Western blot analysis of extracts from LOVO cells, using APBB2 antibody #34412.

Immunohistochemistry analysis of paraffin-embedded human breast carcinoma tissue using APBB2 antibody #34412.



Background

May modulate the internalization of beta-amyloid precursor protein.

The MGC Project Team; Genome Res. 14:2121-2127(2004).

Guenette S.Y., Proc. Natl. Acad. Sci. U.S.A. 93:10832-10837(1996).

Note: This product is for in vitro research use only