

ALS2CR13 Antibody

Catalog No: #34419

Package Size: #34419-1 50ul #34419-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	ALS2CR13 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB
Species Reactivity	Hu
Specificity	The antibody detects endogenous levels of total ALS2CR13 protein.
Immunogen Type	Peptide
Immunogen Description	Synthesized peptide derived from C-terminal of human ALS2CR13.
Target Name	ALS2CR13
Other Names	Amyotrophic lateral sclerosis 2 chromosomal region candidate gene 13 protein; Protein FAM117B; FAM117B;
Accession No.	Swiss-Prot: Q6P1L5NCBI Gene ID: 150864
Uniprot	Q6P1L5
GeneID	150864;
SDS-PAGE MW	60kd
Concentration	1.0mg/ml
Formulation	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

Application Details

Western blotting: 1:500~1:3000

Images



Western blot analysis of extracts from Jurkat cells, using ALS2CR13 antibody #34419.

Background

2 isoforms of the human protein are produced by alternative splicing.

Hillier L.W., Nature 434:724-731(2005).

Hadano S., Nat. Genet. 29:166-173(2001).

Ota T., Nat. Genet. 36:40-45(2004).

Note: This product is for in vitro research use only