

AIG1 Antibody

Catalog No: #34420

Package Size: #34420-1 50ul #34420-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	AIG1 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB
Species Reactivity	Hu
Specificity	The antibody detects endogenous levels of total AIG1 protein.
Immunogen Type	Peptide
Immunogen Description	Synthesized peptide derived from internal of human AIG1.
Target Name	AIG1
Other Names	Androgen-induced protein 1; AIG-1; CGI-103;
Accession No.	Swiss-Prot: Q9NVV5NCBI Gene ID: 51390
Uniprot	Q9NVV5
GeneID	51390;
SDS-PAGE MW	32kd
Concentration	1.0mg/ml
Formulation	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

Application Details

Western blotting: 1:500~1:3000

Images



Western blot analysis of extracts from Jurkat cells and HeLa cells, using AIG1 antibody #34420.

Background

May play a role in androgen-regulated growth of hair follicles.

Seo J., Mol. Cells 11:35-40(2001).

Lai C.-H., Genome Res. 10:703-713(2000).

Ota T., Nat. Genet. 36:40-45(2004).

Note: This product is for in vitro research use only