

ARFIP1 Antibody

Catalog No: #34436

Package Size: #34436-1 50ul #34436-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description


Product Name	ARFIP1 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB IF
Species Reactivity	Hu Ms Rt
Specificity	The antibody detects endogenous levels of total ARFIP1 protein.
Immunogen Type	Peptide
Immunogen Description	Synthesized peptide derived from internal of human ARFIP1.
Target Name	ARFIP1
Other Names	ADP-ribosylation factor interacting 1; arfaptin-1; ARFP1;
Accession No.	Swiss-Prot: P53367NCBI Gene ID: 27236
Uniprot	P53367
GeneID	27236;
SDS-PAGE MW	40kd
Concentration	1.0mg/ml
Formulation	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

Application Details

Western blotting: 1:500~1:3000

Immunofluorescence: 1:100~1:500

Images



Western blot analysis of extracts from Jurkat cells, using ARFIP1 antibody #34436.

Immunofluorescence analysis of A549 cells, using ARFIP1 antibody #34436.



Background

Putative target protein of ADP-ribosylation factor.

Kanoh H., J. Biol. Chem. 272:5421-5429(1997).

Premont R.T., Submitted (JAN-1999) to the EMBL/GenBank/DDBJ databases.

The MGC Project Team; Genome Res. 14:2121-2127(2004).

Note: This product is for in vitro research use only