

## ABCA8 Antibody

Catalog No: #34458

Package Size: #34458-1 50ul #34458-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)


## Description

|                       |                                                                                                                                                    |
|-----------------------|----------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Name          | ABCA8 Antibody                                                                                                                                     |
| Host Species          | Rabbit                                                                                                                                             |
| Clonality             | Polyclonal                                                                                                                                         |
| Purification          | The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.                              |
| Applications          | WB                                                                                                                                                 |
| Species Reactivity    | Hu                                                                                                                                                 |
| Specificity           | The antibody detects endogenous levels of total ABCA8 protein.                                                                                     |
| Immunogen Type        | Peptide                                                                                                                                            |
| Immunogen Description | Synthesized peptide derived from internal of human ABCA8.                                                                                          |
| Target Name           | ABCA8                                                                                                                                              |
| Other Names           | ATP-binding cassette sub-family A member 8;                                                                                                        |
| Accession No.         | Swiss-Prot: O94911NCBI Gene ID: 10351                                                                                                              |
| Uniprot               | O94911                                                                                                                                             |
| GeneID                | 10351;                                                                                                                                             |
| SDS-PAGE MW           | 179kd                                                                                                                                              |
| Concentration         | 1.0mg/ml                                                                                                                                           |
| Formulation           | Rabbit IgG in phosphate buffered saline (without Mg <sup>2+</sup> and Ca <sup>2+</sup> ), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol. |
| Storage               | Store at -20°C                                                                                                                                     |

## Application Details

Western blotting: 1:500~1:3000

## Images



Western blot analysis of extracts from HUVEC cells, using ABCA8 antibody #34458.

## Background

---

ATP-dependent lipophilic drug transporter.

Nagase T., DNA Res. 5:355-364(1998).

Ohara O., Submitted (JAN-2004) to the EMBL/GenBank/DDBJ databases.

Tsuruoka S., Biochem. Biophys. Res. Commun. 298:41-45(2002).

---

Note: This product is for in vitro research use only