

DDX24 Antibody

Catalog No: #34466

Package Size: #34466-1 50ul #34466-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	DDX24 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	IHC
Species Reactivity	Hu
Specificity	The antibody detects endogenous levels of total DDX24 protein.
Immunogen Type	Peptide
Immunogen Description	Synthesized peptide derived from internal of human DDX24.
Target Name	DDX24
Other Names	ATP-dependent RNA helicase DDX24; DD24; DEAD-box protein 24; EC 3.6.1.-;
Accession No.	Swiss-Prot: Q9GZR7NCBI Gene ID: 57062
Uniprot	Q9GZR7
GeneID	57062;
SDS-PAGE MW	97kd
Concentration	1.0mg/ml
Formulation	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

Application Details

Immunohistochemistry: 1:50~1:100

Images



Immunohistochemistry analysis of paraffin-embedded human brain tissue using DDX24 antibody #34466.

Background

ATP-dependent RNA helicase Potential.

Wiemann S., Genome Res. 11:422-435(2001).

Ota T., Nat. Genet. 36:40-45(2004).

Beausoleil S.A., Proc. Natl. Acad. Sci. U.S.A. 101:12130-12135(2004).

Note: This product is for in vitro research use only