

DDX55 Antibody

Catalog No: #34470

Package Size: #34470-1 50ul #34470-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	DDX55 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB
Species Reactivity	Hu Ms
Specificity	The antibody detects endogenous levels of total DDX55 protein.
Immunogen Type	Peptide
Immunogen Description	Synthesized peptide derived from internal of human DDX55.
Target Name	DDX55
Accession No.	Swiss-Prot: Q8NHQ9NCBI Gene ID: 57696
Uniprot	Q8NHQ9
GeneID	57696;
SDS-PAGE MW	68kd
Concentration	1.0mg/ml
Formulation	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

Application Details

Western blotting: 1:500~1:3000

Images



Western blot analysis of extracts from HeLa cells, using DDX55 antibody #34470.

Background

Probable ATP-binding RNA helicase.

The MGC Project Team; Genome Res. 14:2121-2127(2004).

The German cDNA consortium; Submitted (SEP-2004) to the EMBL/GenBank/DDBJ databases.

Nagase T., DNA Res. 7:273-281(2000).

Note: This product is for in vitro research use only