

DHX8 Antibody

Catalog No: #34471

Package Size: #34471-1 50ul #34471-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	DHX8 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB
Species Reactivity	Hu
Specificity	The antibody detects endogenous levels of total DHX8 protein.
Immunogen Type	Peptide
Immunogen Description	Synthesized peptide derived from internal of human DHX8.
Target Name	DHX8
Other Names	DDX8; DHX8; HRH1; PRP22; PRPF22
Accession No.	Swiss-Prot: Q14562NCBI Gene ID: 1659
Uniprot	Q14562
GeneID	1659;
SDS-PAGE MW	150kd
Concentration	1.0mg/ml
Formulation	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

Application Details

Western blotting: 1:500~1:3000

Images



Western blot analysis of extracts from HeLa cells, Jurkat cells and HepG2 cells, using DHX8 antibody #34471.

Background

Facilitates nuclear export of spliced mRNA by releasing the RNA from the spliceosome.

Ono Y., Mol. Cell. Biol. 14:7611-7620(1994).

Ohno M., Genes Dev. 10:997-1007(1996).

Jurica M.S., RNA 8:426-439(2002).

Note: This product is for in vitro research use only