

DQX1 Antibody

Catalog No: #34472

Package Size: #34472-1 50ul #34472-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	DQX1 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB
Species Reactivity	Hu Ms
Specificity	The antibody detects endogenous levels of total DQX1 protein.
Immunogen Type	Peptide
Immunogen Description	Synthesized peptide derived from internal of human DQX1.
Target Name	DQX1
Other Names	ATP-dependent RNA helicase DQX1; DEAQ box 1; DEAQ RNA-dependent ATPase DQX1; EC 3.6.1.-;
Accession No.	Swiss-Prot: Q8TE96NCBI Gene ID: 165545
Uniprot	Q8TE96
GeneID	165545;
SDS-PAGE MW	79kd
Concentration	1.0mg/ml
Formulation	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

Application Details

Western blotting: 1:500~1:3000

Images



Western blot analysis of extracts from COS-7 cells, using DQX1 antibody #34472.

Background

3 isoforms of the human protein are produced by alternative splicing.

Ji W., Mamm. Genome 12:456-461(2001).

Ota T., Nat. Genet. 36:40-45(2004).

Hillier L.W., Nature 434:724-731(2005).

Note: This product is for in vitro research use only