

BAGE2 Antibody

Catalog No: #34475

Package Size: #34475-1 50ul #34475-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	BAGE2 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	IHC
Species Reactivity	Hu
Specificity	The antibody detects endogenous levels of total BAGE2 protein.
Immunogen Type	Peptide
Immunogen Description	Synthesized peptide derived from internal of human BAGE2.
Target Name	BAGE2
Other Names	Cancer/testis antigen 2.2; CT2.2;
Accession No.	Swiss-Prot: Q86Y30NCBI Gene ID: 85319
Uniprot	Q86Y30
GeneID	85319;
SDS-PAGE MW	12kd
Concentration	1.0mg/ml
Formulation	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

Application Details

Immunohistochemistry: 1:50~1:100

Images



Immunohistochemistry analysis of paraffin-embedded human testis tissue using BAGE2 antibody #34475.

Background

Unknown. Candidate gene encoding tumor antigens.

Ruault M., Eur. J. Hum. Genet. 10:833-840(2002).

Ruault M., Genomics 81:391-399(2003).

The MGC Project Team; Genome Res. 14:2121-2127(2004).

Note: This product is for in vitro research use only